





PROKAN PE-6100 Plus



	Result	Unit	Reference
HGB		G/DL	13.5 - 17
RBC		10 ⁶ /uL	4.5 - 6
HCT		Z	40 - 50
MCV		FL	78 - 95
MCH		PG	27 - 32
MCHC		G/DL	32 - 36
RDW-SD		FL	37 - 54
RDW-CV		Z	11.5 - 14.5
PLT		10 ³ /uL	150 - 400
MPV		FL	7.4 - 10.4
PCT		Z	0.10 - 0.20
PDW		Z	10 - 17
WBC		10 ³ /uL	4 - 11

Edit Reference

Auto sampler

- 10.4 LCD Display with the Test items , Reference Value and Three Histograms
- Online Alert for the Abnormal Results
- Touch Screen Operation
- Windows OS
- Easy to Search and Edit Patient Information
- Storage Capacity : 35,000 Sample Results Inclusive of Histograms
- Optional Printout Formats
- Real time Guideline for Easy Operation

Differential	Percentage	ref.Range	Residual value / ul	Ref. Range	
Basophils	0	Z	0-1	0	< 340
Eosinophils	0	Z	1-4	0	< 1000
Neutrophils		Z	40-75	0	5,000 - 21,000
Staf	0	Z	0-5	0	< 1000
Segmented	0	Z	40-70	0	1,500 - 7,500
Lymphocytes		Z	20-45	0	2,000 - 11,500
Monocytes		Z	2-8	0	< 2400

Accurate Results & low sample consumption

- Accurate Analysis Technology With low Sample Consumption
- Venous Mode : 9.6ul
- Capillary Mode : 9.6ul
- Pre-Diluted Mode : 20ul

CBC SoftWare



Exclusive distributor



Add : Zohour District, Zagazig 14512 , Egypt .
 Tel efax : 0552329695
 Mob : 01220108221 - 01024025659 - 01557468221
 E-mail : support@guda-tek.com
 Website : www.guda-tek.com